# RePod SR Packaging



EVERYDAY Clean Beauty. repurchase • refill • recycle

EVERYDAY Clean Beauty.

RePod





rPP Mechanically recycled PCR PP



rPET Chemically Recycled PCR PET

## the DESIGN

The Refillable EBM airless jar, RePod, is an aesthetical design for everyday clean beauty and organic skincare solutions, enabling consumers to repurchase, refill, and recycle the empties after usage.

### the JAR

- ✔ Protect organic ingredients without preservatives.
- ✓ Dispense high-viscosity balm, cream, or hydrating jelly
- ✓ Use 70% plastic less than piston-type airless design.
- ✓ Mono-material, recyclable.

### the LOOK

The ivory-toned overcap, pump, jar & screw cap of repurchases are the resin shade of mechanically recycled PCR, and the crystal clear reusable PET case is chemically recycled PCR.

<u>View decoration options</u> to build a distinguished image.


### EVERYDAY Clean Beauty, RePod

# Consumers Repurchase, Refill, and Recycle.

- 1 Remove the emptied jar.
- 2 Unscrew the cap of the purchased refill.
- 3 Assemble the refill with the reusable pump.
- 4 Toss the emptied jar into recycle bin.



View the catalog



RePod

www.srpackaging.com

New Jersey • Shanghai • Taichung • Tokyo



